Parallel Beam Reconstruction

SS 2022

Y. Huang, F. Wagner, J. Maier, A. Preuhs, C. Syben, A. Maier Pattern Recognition Lab (CS 5)

Friedrich-Alexander University Erlangen-Nuremberg









Topics

Tomography

Projection

Hints for Implementation

Image Reconstruction

Important Methods

Central Slice Theorem

Filtered Backprojection


Hints for Implementation






Basic Principles of Tomography

• $\pi \mathbf{O} \mu \mathbf{O} \sigma$ = tomos = slice










Basic Principles of Tomography (2)

• Idea: Observe object of interest from multiple sides







Topics

Tomography

Projection

Hints for Implementation

Image Reconstruction

Important Methods

Central Slice Theorem


Filtered Backprojection

Hints for Implementation





Projection – Physical Observations




- X-ray Attenuation: $I = I_0 e^{-(\int f(x,y)dI)}$
 - I₀: initial X-ray beam intensity
 - f(x, y): absorption coefficient of material at position (x, y). (x, y) lies on beam line I





Projection – Physical Observations (2)



Observed Signal $I = I_0 e^{-p}$





Projection – Physical Observations (4)




Line Integral Data $p = \log(I/I_0)$






Projection Formation







Projection – Mathematical Formulation



$$p(s,\theta) = \int_{-\infty}^{\infty} \int_{-\infty}^{\infty} f(x,y) \delta(x \cos \theta + y \sin \theta - s) dx dy$$





Topics

Tomography

Projection

Hints for Implementation

Image Reconstruction

Important Methods

Central Slice Theorem


Filtered Backprojection

Hints for Implementation






Parallel Beam Sinogram (Ray_rDriven)








Parallel Beam Sinogram (Ray Driven)







Parallel Beam Sinogram (Ray₋Driven)




Rotation angle and detector index determines one ray.





Parallel Beam Sinogram (Ray₋Driven)




The real detector is equivalent to a detector passing origin.





Parallel Beam Sinogram (Ray_rDriven)




Rotation angle θ and detector index s determines one ray.





Parallel Beam Sinogram (Ray_□Driven)



The ray orientation is orthogonal to the detector orientation, determined by θ . Sampling along the ray for each distance Δt and sum them up as integral.





Topics

Tomography

Projection

Hints for Implementation

Image Reconstruction

Important Methods

Central Slice Theorem

Filtered Backprojection


Hints for Implementation





Reconstruction – Simple Example

Solve the puzzle








Reconstruction – Simple Example

· Solve the puzzle



$$x_1+x_3=7$$

$$x_2 + x_4 = 2$$

$$x_1+x_2=5$$


$$x_3 + x_4 = 4$$





Reconstruction – Simple Example

· Solve the puzzle



$$x_1 + x_3 = 7$$

 $x_2 + x_4 = 2$
 $x_1 + x_2 = 5$
 $x_3 + x_4 = 4$

$$x_1 = 3$$

$$x_2 = 2$$

$$x_3 = 4$$

$$x_4 = 0$$





Reconstruction – Simple Example (2)

· Projection can be formulated in matrix notation

$$\mathbf{P} = \mathbf{AX}$$

$$\mathbf{P} = \begin{pmatrix} 7 \\ 2 \\ 5 \\ 4 \end{pmatrix}, \qquad \mathbf{A} = \begin{pmatrix} 1 & 0 & 1 & 0 \\ 0 & 1 & 0 & 1 \\ 1 & 1 & 0 & 0 \\ 0 & 0 & 1 & 1 \end{pmatrix}, \qquad \mathbf{X} = \begin{pmatrix} x_1 \\ x_2 \\ x_3 \\ x_4 \end{pmatrix}$$





Reconstruction – Simple Example (2)

Solve with matrix inverse?

$$A^{-1}P = X$$


· Common problem size:

$$m{A} \in \mathbb{R}^{512^3 \times 512^2 \times 512}$$
 $512^6 \cdot 4 \; \text{Byte} = 2^{9 \cdot 6} \cdot 2^2 \; \text{B} = 2^6 \cdot 2^{50} \; \text{B}$ $= 64 \; \text{PB} = 65536 \; \text{TB}$





Reconstruction – Example Projection







Reconstruction - Example Backprojection





Image: Zeng, 2009






Reconstruction – Example Backprojection (2)








Reconstruction – Example Backprojection (3)








Reconstruction – Example "Negative Wings"







Reconstruction – Example Reconstruction







Topics

Tomography

Projection

Hints for Implementation

Image Reconstruction

Important Methods

Central Slice Theorem

Filtered Backprojection

Hints for Implementation





Fourier Transform

1D Fourier transform:

$$P(\omega) = \int_{-\infty}^{\infty} p(s) e^{-2\pi i s \omega} \mathrm{d}s$$

• 1D inverse Fourier transform:

$$p(s) = \int_{-\infty}^{\infty} P(\omega) e^{2\pi i s \omega} d\omega$$





Convolution

· Convolution:

$$(f*g)(t) = \int_{-\infty}^{\infty} f(au)g(t- au) \mathrm{d} au = \int_{-\infty}^{\infty} f(t- au)g(au) \mathrm{d} au$$

Convolution theorem:

$$q(s) = f(s) * g(s)$$

$$Q(\omega) = F(\omega) \cdot G(\omega)$$





Topics

Tomography

Projection

Hints for Implementation

Image Reconstruction

Important Methods

Central Slice Theorem

Filtered Backprojection

Hints for Implementation





Central Slice Theorem




Image: Zeng, 2009





Idea for Reconstruction




Image: Zeng, 2009





Topics

Tomography

Projection

Hints for Implementation

Image Reconstruction

Important Methods

Central Slice Theorem

Filtered Backprojection

Hints for Implementation





Filtered Backprojection

- Transform between polar coordinates and Cartesian coordinates, the Jacobian coefficient ω

$$\omega_{x} = \omega \cos \theta, \quad \omega_{y} = \omega \sin \theta$$

Fourier transform in polar coordinates:

$$f(x,y) = \int_0^{2\pi} \int_0^{\infty} F_{\text{polar}}(\omega,\theta) e^{2\pi i \omega (x \cos \theta + y \sin \theta)} \omega d\omega d\theta$$

• Multiplication with $H(\omega) = |\omega|$ in Fourier domain

Convolution with h(s) in image domain





Filtered Backprojection - Practical Algorithm

· Apply Filter on the detector row:

$$q(s,\theta) = h(s) * p(s,\theta)$$

Backproject q(s, θ):

$$f(x,y) = \int_0^{\pi} q(s,\theta)|_{s=x\cos\theta+y\sin\theta} d\theta$$

Backprojection means to add the detector value to each pixel
 → For each pixel at each rotation, you need to find the
 corresponding detector position and read its value.





Discrete Spatial Form of the Ramp Filter

- Find the inverse Fourier transform of $|\omega|$
- Set cut-off frequency of the ramp filter at $\omega = \frac{1}{2}$





Discrete Spatial Form of the Ramp Filter

- Find the inverse Fourier transform of $|\omega|$
- Set cut-off frequency of the ramp filter at $\omega=rac{1}{2}$

$$h(s) = \int_{-\frac{1}{2}}^{\frac{1}{2}} |\omega| e^{2\pi i \omega s} d\omega$$





Discrete Spatial Form of the Ramp Filter

- Find the inverse Fourier transform of $|\omega|$
- Set cut-off frequency of the ramp filter at $\omega = \frac{1}{2}$

$$h(s) = \int_{-\frac{1}{2}}^{\frac{1}{2}} |\omega| e^{2\pi i \omega s} d\omega$$

$$h(s) = \frac{1}{2} \frac{\sin \pi s}{\pi s} - \frac{1}{4} \left[\frac{\sin \left(\frac{\pi s}{2} \right)}{\frac{\pi s}{2}} \right]^2$$





Discrete Spatial Form of the Ramp Filter (2)

$$h(s) = \frac{1}{2} \frac{\sin \pi s}{\pi s} - \frac{1}{4} \left[\frac{\sin \left(\frac{\pi s}{2} \right)}{\frac{\pi s}{2}} \right]^2$$





Discrete Spatial Form of the Ramp Filter (2)

$$h(s) = \frac{1}{2} \frac{\sin \pi s}{\pi s} - \frac{1}{4} \left[\frac{\sin \left(\frac{\pi s}{2} \right)}{\frac{\pi s}{2}} \right]^2$$

• Convert to discrete form: Let $s = n \cdot s_0$ (*n* integer, s_0 spacing)

$$h(n \cdot s_0) = \left\{ egin{array}{ll} rac{1}{4s_0^2} & n = 0 \ 0 & n ext{ even} \ -rac{1}{n^2\pi^2s_0^2} & n ext{ odd} \end{array}
ight.$$





Discrete Spatial Form of the Ramp Filter (2)

$$h(s) = \frac{1}{2} \frac{\sin \pi s}{\pi s} - \frac{1}{4} \left[\frac{\sin \left(\frac{\pi s}{2} \right)}{\frac{\pi s}{2}} \right]^2$$

• Convert to discrete form: Let $s = n \cdot s_0$ (*n* integer, s_0 spacing)


$$h(n \cdot s_0) = \left\{ egin{array}{ll} rac{1}{4s_0^2} & n = 0 \ 0 & n ext{ even} \ -rac{1}{n^2\pi^2s_0^2} & n ext{ odd} \end{array}
ight.$$

Also known as the "Ramachandran-Lakshminarayanan" convolver





Discrete Spatial Form of the Ramp Filter (3)







Topics

Tomography

Projection

Hints for Implementation

Image Reconstruction

Important Methods


Central Slice Theorem

Filtered Backprojection

Hints for Implementation















Get each pixel position x, y.








According to rotation angle θ , determine detector orientation.








Want to calculate the position s, where the ray hits the detector.







s can be calculated simply as the vector p projected to the direction vector v.